



## Product Description

The KAPALongRange PCR system is a blend of *Taq* DNA polymerase and a modified archaeal (Type B) DNA polymerase possessing proofreading capability. This two-enzyme system is designed specifically to support long range and/or sensitive PCR. The KAPALongRange system polymerizes DNA from a primer annealed to a DNA template in the presence of deoxyribonucleotide triphosphates. Both enzymes possess 5'→3' polymerase activity, but only *Taq* possesses double strand dependent 5'→3' exonuclease activity and only the Type B polymerase possesses 3'→5' exonuclease (proofreading) activity.

KAPALongRange possesses higher fidelity than *Taq* polymerase.

The majority of PCR products from a KAPALongRange DNA Polymerase reaction have a single deoxyadenosine added to the 3' ends of the amplicon and can be used for ligation into single 3'-T-overhang cloning sites.

All KAPALongRange kits are supplied with a 5x reaction buffer. This buffer has been formulated to support long range and low template applications. Magnesium chloride solution is supplied separately to accommodate PCR optimization.

## Product Applications

The KAPALongRange PCR system is ideally suited for:

- PCR amplification of long targets and/or PCR using low amounts of template DNA
- Standard short- and mid-range PCR amplification
- Production of PCR products to be used for ligation into 3'-T-overhang cloning vectors

## Product Performance


### LONG RANGE AND HIGH SENSITIVITY

KAPALongRange is able to amplify long templates from limiting complex genomic DNA.

**Figure 1:** Amplification of 15 kb fragment (tPA gene) from human chromosomal DNA. Lanes: (M) Marker, (1) 25 ng human gDNA, (2) 5 ng human gDNA, (3) 1 ng human gDNA. 35 cycles, 68°C extension temperature, 2.5 units of enzyme.

KAPALongRange exhibits high yield, greater sensitivity and specificity compared with leading competitive long range kits.

**Figure 2:** Amplification of 15 kb fragment from 10 pg, 10 fg, and 1 fg of lambda DNA. Lanes: (1-3) KAPALongRange, (M) Marker, (4-6) Competitor R, (7-9) Competitor T.



next generation thinking  
in enzyme technology




Long Range. High Sensitivity.

# >> KAPALongRange | DNA Polymerase

## Product Performance

KAPALongRange can replace *Taq* polymerase using standard PCR cycling conditions and less enzyme per reaction. KAPALongRange amplifies short-range and mid-range targets with greater sensitivity and specificity as compared with *Taq* polymerase.

Figure 3



**Figure 3:** Amplification of 6 log dilution series of lambda DNA starting with 1 ng to 10 fg. Amplicons were 2 kb, 8 kb, and 15 kb in length. 35 cycles, 72°C extension temperature, 1.25 units of enzyme.

## KAPALongRange Kit Components:

- KAPALongRange DNA Polymerase (5 U/μl in storage buffer)
- 5x KAPALongRange Reaction Buffer
- 25 mM MgCl<sub>2</sub> solution
- 250μl KAPA dNTP Mix (10mM each)

KAPALongRange	
Product Code	Kit Size
KK3003	250 units
KK3004	100 units

For custom orders please contact: [sales@kapabiosystems.com](mailto:sales@kapabiosystems.com)

### Boston, Massachusetts, United States

600 West Cumming Park, Suite 5350 Woburn, MA, 01801 U.S.A.  
Tel: +1 781 497 2933 Fax: +1 781 497 2934  
Email: [info@kapabiosystems.com](mailto:info@kapabiosystems.com)

### Cape Town, South Africa

Research, Development, and Manufacturing Facility  
2nd Floor, Old Warehouse Building, Black River Park, Fir Road,  
Observatory, 7925 Cape Town, South Africa  
Tel: +27 21 448 8200 Fax: +27 21 448 6503  
Email: [info@kapabiosystems.com](mailto:info@kapabiosystems.com)

