

**Stain your nuclei  
in 5 minutes...**

**abcam**<sup>®</sup>  
*discover more*

**DRAQ5**<sup>™</sup>

## **Rapid DNA fluorescent dye**

- Rapid staining of dsDNA/nuclei of LIVE and fixed cells
- Easy to use – no lyse, no wash, no RNase needed
- Ideal for use with GFP, FITC labels and DyLight<sup>®</sup> 488
- $Ex\lambda_{max}$  646 nm,  $Em\lambda_{max}$  681 nm / 697 nm intercalated with dsDNA
- No photobleaching effect
- Fully compatible with CyGEL<sup>™</sup> thermoreversible gel (ab109204)



Discover more at [abcam.com/DRAQ5](http://abcam.com/DRAQ5)